



RespSense™




Capnograph

MONITORING THE PATIENT



1. Visually inspect the monitor and sample line to make sure there are no visible signs of damage.
2. Use a new sample line, Nafion® tubing (when used), moisture trap and filter for each patient.
3. Connect the sample line to the Luer lock connector on the moisture trap.
4. Apply the sample line to the patient.
5. Turn on the monitor by pressing the ON/OFF  button until you hear a beep.
6. Verify all parameters are displayed.
7. Turn off the monitor by pressing the ON/OFF  button for more than 1 second.



ALARMS


1. At start-up, alarms are automatically silenced for approximately 2 minutes.
2. Alarm limits are automatically set to the factory defaults at each start-up.
3. **Change Alarm Settings** – Increase alarm limits with ; decrease with .
4.  **Audible Alarm Pause/Resume button** – Pause/activate the audible alarm. If paused, it will remain silent for approximately 2 minutes unless the button is pressed again or a new audible alarm occurs.

POWER AND BATTERY

1. Plug the LifeSense power supply into the monitor or operate on battery power.
2. The green charge indicator  comes on as soon as the monitor is connected to the power outlet. It indicates the battery is charging.
3.  Indicates remaining battery capacity.
4. When fully charged, the LifeSense can run on battery for approximately 8 hours.

REPLACING THE MOISTURE TRAP AND FILTER

1. Place the filter in the moisture trap as shown (1). Be sure to orient the filter properly inside the opening at the top of the moisture trap with the hydrophobic material toward the trap.
2. Slide the moisture trap into position (2) using the guide marks on the back of the monitor. Align the groove at the top of the trap with the ridge at the top of the trap housing.
3. Press the moisture trap into position using the tab (3). The trap is placed in top first, then bottom. There will not be an audible or tactile click but the trap should be flush with the monitor housing.
4. To remove the moisture trap and replace the filter, reverse the three steps above. Pull out trap using the bottom tab and remove filter (may still be inside of trap housing).



THIS GUIDE IS FOR REFERENCE ONLY. Refer to the Operator's Manual for complete information, warnings, cautions and specifications for the product you are using.

Visit nonin.com for more information and video tutorials.

R_x only



nonin.com

Nonin Medical, Inc.
13700 1st Avenue North
Plymouth, MN 55441-5443 • USA
Tel: +1.763.553.9968
+1.800.356.8874
Fax: +1.763.577.5521
E-mail: info@nonin.com

Nonin Medical B.V.
Prins Hendriklaan 26
1075 BD Amsterdam • Netherlands
Tel: +31 (0)13-79 99 040
Fax: +31 (0)13-79 99 042
E-mail: infointl@nonin.com