



9847V

EtCO₂ Detector/ Digital Pulse Oximeter

High-Tech Simplicity designed for Veterinary Use

Ideal for surgery or dental procedures on intubated patients, the unique 9847V combines semi-quantitative end-tidal carbon dioxide detection with digital pulse oximetry into one portable handheld monitor for maximum functionality. This vital tool is useful to aid and verify endotracheal tube placement throughout the intubation period. During the procedure, an audible breath tone is present. Additionally, the 9847V provides accurate blood oxygen saturation, pulse rate data, and a visual respiration indicator. Truly versatile, the 9847V allows simultaneous or independent use of the EtCO₂ detector and pulse oximeter.

Product Highlights

- ✓ Easy-to-Use - Bright LED displays and intuitive icons
- ✓ "First Breath" Active - No warm-up required
- ✓ Cost Effective - No calibration needed and long battery life
- ✓ Breath Tone - Audible indication of detected breath
- ✓ Indication of Pulse Quality - Tricolor perfusion indicator
- ✓ Veterinary Specific - Hardware and software
- ✓ Durable - Rugged construction to withstand daily use
- ✓ Compact Size - Easily portable; requires little space for storage

✓ Where People Make the Difference



Pure Functionality

The powerful 9847V provides “first breath” activation by eliminating any warm-up time – simply turn on the power and start taking measurements. NONIN’s patented semi-quantitative EtCO₂ technology, combined with mainstream sampling, eliminates the need for calibration, thus reducing operating expenses.



The sensor includes a fully patented airway adapter tube with a low 6 cc dead space and anti-fog coating for superior performance. An easy-to-interpret and highly visible LED bar graph is standard on all models.

Combined with NONIN’s PureLight® sensors, the 9847V delivers the outstanding performance you need to get the job done in a demanding environment.



2000SL – Clip on sensor for tongue applications, paw pads and well-vascularized areas.



2000T – For placement on the underside base of the tail or other well vascularized surfaces.



2000SA – Wrap sensor for placement on the toe (large animal) or base of the tail (small animal).

Features & Alarm Functions


9847V
Full-Alarm



Features

| | |
|----------------------------------|---|
| Breath Tone | ✓ |
| Adjustable Volume | ✓ |
| Variable Pitch - CO ₂ | ✓ |
| Alarm Functions | |
| Audible Mute | ✓ |
| Audible Disable | ✓ |
| Apnea | ✓ |
| SpO ₂ High/Low | ✓ |
| Pulse Rate High/Low | ✓ |

Please contact your authorized distributor for more information.



© 2005 Nonin Medical, Inc.
Nonin®, the Nonin Logo® and PureLight® are registered trademarks of Nonin Medical, Inc.
P/N 2937-000-03

Specifications

Oximeter

| | |
|--|---|
| Oxygen Saturation Range (%SpO ₂) | 0% to 100% |
| Pulse Rate Range | 18 to 450 pulses per minute |
| CO ₂ Range | 0 to >75 mmHg |
| CO ₂ Response Time | 250ms |
| Respiration Rate Range | 1 to 60 breaths/minute |
| Breath Detection Threshold | 5 mm Hg |
| Measurement Wavelengths | 660 nanometers 910 nanometers 4260 nanometers |
| Displays | Two 3-digit 7-segment LEDs (light emitting diodes) Tricolor LED segments 8 segment LED bar graph |
| Accuracy | |
| Blood Oxygen Saturation (%SpO ₂ ±1 SD)* | 70 - 100% ± 3 digits for the Model 2000SL, 2000SA and 2000T sensors, below 70% is not specified for all sensors |
| Pulse Rate | ±3% ±1 digit |
| Temperature | |
| Operating (pulse oximeter) | -4° to +122°F (-20° to +50°C) |
| Operating (CO ₂ detector) | +32° to +122°F (0° to +50°C) |
| Storage | -22° to +122°F (-30° to +50°C) |
| Humidity | |
| Operating | 15% to 90% noncondensing |
| Storage | 10% to 95% noncondensing |
| Power Requirements | 6AA alkaline batteries |
| Battery Life | |
| Operating | 90 hours typical operation (pulse oximeter alone) 20 hours typical operation (CO ₂ detector and pulse oximeter) |
| Dimensions | 3 x 6 x 1 in. (8 x 15 x 2 cm) |
| Weight | 11 oz. with batteries (310 g) |
| Memory | 24 hours (SpO ₂ , PR, breath interval) |
| Warranty | 3 years |

Classifications per IEC 60601-1/UL60601-1

| | |
|----------------------|---------------------------------------|
| Type of Protection | Internally powered (on battery power) |
| Degree of Protection | Type BF-Applied Part |
| Mode of Operation | Continuous |

* S.D. (Standard Deviation) is a statistical measure; up to 32% of the readings may fall outside of these limits.

Quality Systems are registered to ISO 13485:2003 and EC certificates per Annex II, clause 3 of EC Directive No. 93/42/EEC concerning medical devices.



Specifications are subject to change without notice.

Your Partner of Choice



With over 20 years of experience and dedication to the design and support of noninvasive monitoring devices, NONIN has helped many medical professionals meet clinical and economic objectives. Trusted by clinicians worldwide, countless tests are performed every day on thousands of NONIN pulse oximeters in more than 125 countries. NONIN's dedication to technological leadership, precision manufacturing and uncompromised customer support ensures quality products and service you expect.



Nonin Medical, Inc.
13700 1st Avenue North
Plymouth, MN 55441-5443, USA
Phone: +1.763.553.9968 Toll Free: 1.800.356.8874
Fax: +1.763.553.7807
info@nonin.com

✓ Where People Make the Difference