



Nonin Medical, Inc.  
13700 1st Avenue North  
Plymouth, MN 55441-5443  
USA

## MR COMPATIBILITY CERTIFICATION

Nonin Medical Inc certifies that Model 7500FO Digital Pulse Oximeter can operate in a magnetic resonance (MR) environment with a 3T (Tesla) magnet, outside the 200 Gauss line.

Testing in support of this certification is on file at Nonin Medical, in the form of report number QATR5602-000.

*Kim E. Aves*

Signature

*03 June 2005*

Date

Kim E. Aves, Regulatory Affairs Specialist

Printed Name and Title