


Introducing Upgrades to WristOx₂[®] Bluetooth[®] Technology

WristOx₂ Model 3150 with Bluetooth Low Energy (BLE)
Set to Replace Original Model 3150



9 November 2018

Bluetooth Classic is being discontinued. Here’s what that means for the WristOx₂ and for you.

The original WristOx₂[®] Pulse Oximeter Model 3150 uses a Bluetooth transceiver that will not be made after March 2019. To continue to supply wireless pulse oximetry solutions, Nonin is introducing a new generation in the WristOx₂ line: Model 3150 with Bluetooth Low Energy (BLE). We recommend transitioning to the Model 3150 with BLE for your wireless pulse oximetry needs.

The Model 3150 with BLE incorporates a low energy Bluetooth transmitter, both improving wireless connectivity and extending the battery life of the device. Other key updates include:

- Improved connectivity:** Memory download via USB or Bluetooth LE makes the upgraded device convenient in most patient care settings.
- Longer battery life:** Two standard AAA batteries provide a minimum of 44 hours of use.
- Free Mobile App:** The NoninConnect application*, available through the Apple store, allows users to track SpO₂ readings and share them by email.
- Affordable Price:** The new Model 3150 with BLE will be available for the same price as the original Model 3150.

WristOx₂ Product Line

Nonin’s updated WristOx₂ product line will continue to include the Model 3150 with USB in addition to the newly introduced Model 3150 with BLE. The Model 3150 with USB is recommended for situations where wired connections are preferred and is offered at a lower price points compared to the BLE version. To distinguish the devices, we are introducing a blue colored lens to the Bluetooth device.

	3150 USB	3150 with BLE
Bluetooth	No	4.2
Battery Life	53 hours	44 hours
Bluetooth Range (m)	N/A	60
nVision Software Compatible	Yes (USB cable only)	Yes (USB cable only)
NoninConnect App (iOS)	No	Yes (Bluetooth)



Model 3150 with USB Model 3150 with BLE

Ordering Information

Contact Nonin’s customer service or your account manager to place an order.

SKU	Description	SKU	Description
3150BLE-0101	3150 with BLE, 8000SM-WO2	3150USB-0101	3150 USB, 8000SM-WO2
3150BLE-0102	3150 with BLE, 8000AA-WO2	3150USB-0102	3150 USB, 8000AA-WO2
3150BLE-0103	3150 with BLE, 8000J-WO2	3150USB-0103	3150 USB, 8000J-WO2
3150BLE-0104	3150 with BLE, 8000SS-WO2	3150USB-0104	3150 USB, 8000SS-WO2
3150BLE-0105	3150 with BLE, 8000SL-WO2	3150USB-0105	3150 USB, 8000SL-WO2

Nonin will continue to support any original WristOx₂ 3150 devices currently in use. We will contact you again before the original WristOx₂ 3150 is discontinued—and we encourage you to switch to the new WristOx₂ 3150 with BLE or WristOx₂ 3150 USB with your next order.

* The NoninConnect app is not cleared as a medical device

To place an order, contact Nonin Customer Service

Nonin Medical, Inc.
13700 1st Avenue North
Plymouth, MN 55441-5443 USA
Tel: +1.763.553.9968
800.356.8874
Fax: +1.763.577.5521
E-mail: info@nonin.com

Nonin Medical B.V.
Doctor Paul Janssenweg 150
5026 RH Tilburg
The Netherlands
Tel: +31 (0)13 - 79 99 040
Fax: +31 (0)13 - 79 99 042
E-mail: infointl@nonin.com